Practicon® E-Prop™ Mouth Prop or Soft-Prop Mouth Prop

71-02411, 71-02412, 71-02413

Instructions

1. Remove from packaging and sterilize according to separate instructions before use.

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2. When ready to use, insert suction tube into piece (Figure 1).

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3. Pull until suction tube locks in place (Figure 2).

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4. Place in mouth over back molars (Figure 3).

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5. Do not place over front teeth (Figure 4).

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6. After use, rinse and inspect for damage. If damaged, discard. If not damaged, it can be sterilized and used again.

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Figure 3


Figure 4

Autoclaving Procedures for the E-PROP

The E-Prop is reusable and must be sterilized after each use. Autoclaving is the most dependable and preferred method of decontamination and sterilizing of the E-Prop.

All items should be carefully cleaned before autoclaving to prevent baking contaminants onto the surface of the plastic. After cleaning, the E-Prop should be rinsed thoroughly in distilled water before autoclaving. Certain chemicals which have no appreciable effect on plastics at room temperature may cause deterioration at autoclaving temperatures and therefore must be removed.

Preparation and loading of materials

Use a biohazard or autoclavable bag. These bags are available in a variety of sizes and some are printed with an indicator that changes color when sterilized. Place the E-Prop in an appropriate sized biohazard or autoclave bag. Place the bags in a rigid autoclavable container during autoclaving. Follow the bag manufacturer's instructions. Recommended autoclave cycle is 121°C, 15 psig (1 bar) for 20 minutes. Use fast exhaust and dry cycle.

Chemical Sterilants

Sterilization implies the killing of all living organisms. Chemical sterilants can be used to decontaminate the E-Prop. It is important to choose a chemical sterilant that has been approved and proven effective for this use. Follow manufacturer's instructions. Some chemical sterilants may discolor the E-Prop. Therefore, we recommend autoclaving as the preferred method of sterilization.

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Manufactured for Practicon, Inc. by: Better Dental Products, 8540 "I" Street, Omaha, Nebraska 68127